

Personalized Molecular Medicine™

## Complete Suite of NGS Clonality Assays with Bioinformatics - Identification and Tracking Patient-Specific Clones


Presha Shah, Ph.D.  
Development Scientist

## Outline


- Part I: LymphoTrack Products
  - What are different types clonality targets?
  - How diversity is generated?
  - What is LymphoTrack Assay used for?
  - Why test clonality?
  - How clonality was tested?
  - Overview of *TRB* gene, Multiplexing and Workflow
  - LymphoTrack Assay for MRD
- Part II: LymphoTrack Data Analysis – Identification and Tracking Patient- specific Clones

## Clonality Targets


### B-Cell Receptors (BCR)/Immunoglobulins (Ig)




### T-Cell Receptor (TCR)




## Human *IGH* locus




## Human *IGK* locus



## Human *TRB* locus



## Human *TRG* locus



Ginsburg, C. et al. (2011, Feb 8). *IMGT Repertoire (IG and TR)*. Retrieved from <http://imgt.org/>


invivoscribe

7

## Generation of Diversity: *IGH* Rearrangement

### Gene Rearrangements Are

- Lymphoid-specific
- Developmentally-regulated
- Sequential



*IGH*, *TRB* & *TRD*: V, D, and J  
*IGK*, *IGL*, and *TRG*: V and J regions

invivoscribe

8

## B- and T-Cell Gene Rearrangements

- During B- and T-Cell Development and Maturation
- Genetic Recombination Process at the DNA Level
- Each of the V, D and J Gene Segments are Randomly Recombined
- Encode Unique Antigen Receptors → up to  $10^{14}$
- Sources of Diversity:
  - Gene Rearrangements
  - N-region Diversity
  - Somatic Hypermutations

invivoscribe

9

## Polyclonal vs. Clonal Progression

### Polyclonal Progression



### Clonal Progression



Highly indicative of B- or T-cell malignancy

invivoscribe

10

## Why Test for B- and T-Cell Clonality?


Leukemias and lymphomas can be challenging to diagnose by

- morphology
- immunohistochemistry
- flow cytometry

Using these methods, 5-15% of cases result in inconclusive diagnoses.

Diagnosis of lymphoid malignancies is greatly supported and facilitated by clonality testing.

- Being adopted in routine diagnostics specifically to track disease with MRD testing.



invivoscribe

11

## Our Focus

PCR-based molecular testing for B- and T-cell malignancies associated with leukemias and lymphomas

- Gene Clonality / Rearrangements
- Chromosome Translocations
- Gene Mutations

RUO (Research Use Only)


CE Marked IVD (Clinical Diagnostic Utility)

GPR controls (General Purpose Reagents (DNA and RNA))


invivoscribe

12


## Detection Methods



» Gel  
IdentiClone™ Assays



» ABI  
IdentiClone™ Assays



» NGS  
LymphoTrack® (Dx) Assays

invivoscribe

## Evolution of Clonality Testing

With the evolution of technology, clonality analysis can now be performed using Next Generation Sequencing (NGS) technologies

PCR-based clonality assays are now widely accepted as the gold standard, with improved detection limits

Fragment Analysis

NGS

Limit of Detection

- Sample type suitability

invivoscribe

## Advantages of NGS Clonality Testing

- Determines the **DNA sequence** of clonal rearrangements
- Reduces **subjectivity in interpretation**
- Identifies the **full range** of clonal populations in a specimen
- Allows both **identification and tracking** of clonal populations with the **same reagents and workflow**
- Can be **standardized between labs** and **cleared or approved by Regulatory Agencies Worldwide**

invivoscribe


## LymphoTrack® Dx Portfolio

### Available Assays & Software

- TRG
- TRB - **New!** MiSeq
- TRB - (Ion PGM/SS in development)
- IGHV Leader
- IGH FR1/2/3 Combo
- IGH FR1
- IGH FR2
- IGH FR3
- IGK
- LymphoTrack® Dx Software
- LymphoTrack® MRD Software - Research Use Only

invivoscribe

## LymphoTrack® & LymphoTrack® Dx Workflow




DNA → Amplification → Purify PCR Products → Library Quantification → Emulsion PCR & Enrichment → Run Prep → FASTQ files → Analyze on PC (3-5 min/sample)

Run Prep: 300/500 cycle v2 kit, 600 cycle v3 kit

Run Prep: 24/39 h, 55 h

invivoscribe

## One assay – Multiple applications!



- Clonality
- Somatic Hypermutation
- Minimal Residual Disease\*
- CAR-T Cell & Immunotherapy Tracking

invivoscribe

## LymphoTrack® & LymphoTrack®Dx Advantages

- **Tests Compatible with Range of Specimen Types**
- **ONE-STEP PCR Master Mixes**
- **MULTIPLEXING – Reduces Costs**
  - Combine Multiple Samples - up to 12 (Ion PGM) & 24 (MiSeq)
  - Combine Multiple Assays - up to 72 or 168 samples
- **Unparalleled Sensitivity**
  - Unparalleled clonality detection with the ability to identify and track the specific sequence of clonal populations for MRD studies.
- **INCLUDED Analysis SOFTWARE Package**

## LymphoTrack®Dx *TRB* Assay- MiSeq®



## LymphoTrack®Dx *TRB* Assay- MiSeq®




- ❑ T-cell malignancies arise from clonal expansion of single cell.
- ❑ During early T-cell development, somatic rearrangements occur within T cell receptor beta (*TRB*) locus that bring together, sequentially, the joining of D and J gene segments, followed by joining of a V segment to the DJ pair
- ❑ Next-generation sequencing (NGS) based gene rearrangement methods may improve the sensitivity of clonal detection and identify the specific V(D)-J DNA sequences
- ❑ T cell receptor gamma (*TRG*) locus rearranges prior to T cell receptor beta (*TRB*) locus and testing both loci can identify the vast majority of T-cell malignancies.


## IdentiClone® *TCRB* Gene Clonality Assays Capillary Electrophoresis (Size)

• Tube A: V + J (240-285 nt)

• Tube B: V + J (240-285 nt)



### 1 step PCR



- ❑ 24 indices for Panel and 8 indices for Kit A
- ❑ Positive and Negative control

Kit	Name	# of Indices/Kit	Indices	# of Runs	# of Controls	# of Boxes
Panel	LymphoTrack® Dx TRB Assay Panel – MiSeq®	24	1 - 27	5	3 + 3	1
Kit A	LymphoTrack® Dx TRB Assay Kit A – MiSeq®	8	1 - 8	5	3 + 3	1

## Multiplex- One Step PCR

### Primers

Vβ + Db + Jβ

### DNA

- Peripheral Blood (PB)
- Formalin Fixed Paraffin Embedded (FFPE)
- Bone Marrow (BM)

### Master Mix

- Master Mix
- dNTP
- Polymerase

### PCR

Upto 24 Indices


Pooling & Sequencing

MiSeq

Sequence up to 24 Samples per Target

Multiplex different target with same index


## TRB Positive (Clonal)



invivoscribe

25

## TRB Negative (Non-Clonal)



invivoscribe


26

## Minimal Residual Disease (MRD)

Leukemic cells that remain during or after treatment when a patient is in remission

### Major cause of relapse; Use to determine:

- Has treatment eradicated clonal cells?
- Efficacy of different treatments
- Monitor patient remission status
- Detect recurrence



### Benefits of NGS vs. Flow Cytometry & ASO-PCR:

- Better sensitivity & specificity
- Can be standardized across treatment centers
- No custom assay development required
- Sequence identity / frequency distribution of all alleles
- Track multiple clones

invivoscribe

27

## LymphoTrack® Assays & MRD

- Same LymphoTrack Assay Kits As *IGH*, *IGK*, *TRG*, *TRB* Clonality
- Same Library Preparation
- Sensitivity Only Limited By DNA Input In PCR Amplification Step



invivoscribe

MRD Applications are For Research Use Only. Not for use in Diagnostic Procedures.

## Summary of All LymphoTrack Assays - MiSeq


	MiSeq Primers					
	IGHV Leader 5'NM	IGHV FR1	IGHV FR2	IGHV FR3	IGK	TRG
Target Size (bp)	483	295	243	358	222	147
Amplification Size Inc. Target, Index, and Adapters (bp)	650	450	390	550	410	300
Recommended Sequencing Kit*	MiSeq v1 Reagent (300-cycle)	MiSeq v2 Reagent (300-cycle) or MiSeq v1 Reagent (300-cycle)	MiSeq v2 Reagent (300-cycle) or MiSeq v1 Reagent (300-cycle)	MiSeq v2 Reagent (300-cycle) or MiSeq v1 Reagent (300-cycle)	MiSeq v2 Reagent (300-cycle) or MiSeq v1 Reagent (300-cycle)	MiSeq v2 Reagent (300-cycle) or MiSeq v1 Reagent (300-cycle)

When multiplexing amplicons of different gene targets it is important to use the appropriate sequencing chemistry. The number of sequencing cycles must be sufficient to sequence the largest amplicon in the multiplex.

invivoscribe

29

## Multiplexing Can Reduce Run Cost




invivoscribe

30

**invivoscribe**  
Personalized Molecular Medicine™

**LymphoTrack® Data Analysis – Identification and Tracking Patient-Specific Clones**


Kasey Hutt, Ph.D.  
Staff Scientist,  
Bioinformatics



**LymphoTrack® (Dx) Software Features**

- Simultaneously analyze multiple targets
- Rich data
- Free software upgrades
- Clear, exportable data/graphics
- For MiSeq™, PGM™, S5 .FASTQ files
- Free of charge
- Easy to run
- Local analysis/No internet
- Low computer requirements

**invivoscribe**




**LymphoTrack® Dx Software – CD Content**

**invivoscribe**


**LymphoTrack® Dx Software – Processing Data**

**invivoscribe**

## LymphoTrack® Dx Software – License Agreement




## LymphoTrack® Dx Software – Processing Data




Select target(s) represented in pooled library

## Data Output – TRB Assay



## Importance of Merge Read Summaries

Amplification & Sequencing May Cause Artificial Nucleotide Changes



Merged Read Summary:

The Read Summary tab only shows the top 10 reads after merging with the top 500 reads that differ in 1 or 2 nucleotides

## Where to Find SHM Data?

Sample	SI	471,441								
Rank	Sequence	Length	Raw count	V-gene	J-gene	% total reads	Correlation	Mutation rate to parent (V-gene %)	In-frame (Y/N)	No Stop codon (Y/N)
1	OTCTGTGATGAC	255	4245	JH94-46.01	JH94-52	0.26	0.24	0.00	Y	Y
2	OTCTGTGATGAC	255	542	JH94-46.01	JH94-52	0.08	0.32	0.44	Y	Y
3	OTCTGTGATGAC	255	535	JH94-46.01	JH94-52	0.08	0.40	0.44	Y	Y
4	OTCTGTGATGAC	255	537	JH94-46.01	JH94-52	0.08	0.47	0.44	Y	Y
5	OTCTGTGATGAC	255	490	JH94-46.01	JH94-52	0.07	0.54	0.44	Y	Y
6	OTCTGTGATGAC	255	449	JH94-46.01	JH94-52	0.07	0.61	0.00	Y	Y
7	OTCTGTGATGAC	255	448	JH94-46.01	JH94-52	0.07	0.68	0.00	Y	Y
8	OTCTGTGATGAC	255	394	JH94-46.01	JH94-52	0.06	0.74	0.44	Y	Y
9	OTCTGTGATGAC	255	271	JH94-46.01	JH94-52	0.04	0.78	0.44	Y	Y
10	OTCTGTGATGAC	255	252	JH94-46.01	JH94-52	0.04	0.82	0.44	Y	Y
11	OTCTGTGATGAC	255	236	JH94-46.01	JH94-52	0.03	0.85	0.44	Y	Y
12	OTCTGTGATGAC	255	199	JH94-46.01	JH94-52	0.03	0.88	0.44	Y	Y
13	GGCTTGGATGAC	300	131	JH94-35.01	JH94-52	0.03	0.89	0.56	Y	Y

## LymphoTrack® Dx Visualization Automation

PDF Report Automatically Generated by Analysis Software for Each Target and Sample

Report will include

- Top 10 merged read summary
- Top 200 sequence frequency graph
- Top 200 gene usage graph
- Top 200 read summary





## MRD Software - Data Input

**Step 5: Entering Replicates**

5.1 Select Unique Reads File  
5.2 Entering Replicate Information: DNA Amount, Number of Replicates and Save Replicates

Replicate	Amount of DNA	Unique Reads File	Number of Reads	Edit	Delete
1	700	H:\Innovoscribe\Lab\...	27036	[Edit]	[Delete]

invivoscribe

MRD Software is for Research Use Only (RUO)

49

## MRD Software - Data Input

**Step 6: Save Replicates**

6.1 Click "Back"

Replicate	Amount of DNA	Unique Reads File	Number of Reads	Edit	Delete
1	700	H:\Innovoscribe\Lab\...	27036	[Edit]	[Delete]

invivoscribe

MRD Software is for Research Use Only (RUO)

50

## MRD Software - Data Input

**Perform Analysis**

Project Name: Test1  
Sequence Name: SeqBank1

invivoscribe

MRD Software is for Research Use Only (RUO)

51

## MRD Software - Detected Sequence

**Detected Sequence 1**

Sequence name: SeqBank1  
Sequence input: CATTGTGATACACTTCAACAGCTATATATGACTGCTGACAGGCCCCCTGACAGAGGCTTGAATGATGGAAAT...

invivoscribe

MRD Software is for Research Use Only (RUO)

52

## MRD Software - Not Detected Sequence

**Not Detected Sequence 1**

Sequence Not Detected checked: 1 Replicates  
Sequence input: CATTGTGATACACTTCAACAGCTATATATGACTGCTGACAGGCCCCCTGACAGAGGCTTGAATGATGGAAAT...

invivoscribe

MRD Software is for Research Use Only (RUO)

53

## Project Planner

**Project Planner**

# of Replicates: 2  
# of Reads per Replicate: 27036  
Total Reads per Replicate: 54072

invivoscribe

MRD Software is for Research Use Only (RUO)

## PDF Output

## Future Additions

- Better quantification
- Optimized data flow between Clonality and MRD analysis
- Timepoint aggregation

## Thank You Everyone at IVS!

### Assay Development Team

- Adam Bouzaflour
- Austin Jacobsen
- Brandon Glivens
- Edgar Vigil
- Jeff Panganiban
- Maggie Kaminsky
- Presha Shah
- Ying Huang
- Wenhui Huang
- Martin Blankford
- Selena Zheng
- Jean Xie
- Jeff Miller
- Sam An

### Bioinformatics/SW Dev Team

- Kasey Hunt
- Jamie Lantry
- Joshua Waldman
- Nathan Nichols-Roy

### Marketing/Product Development Team

- Kevin Dobyns
- Zephry Fure
- Rosanne Roncarolo de Vries
- Oscar Rodriguez, Jr.

LabPMM (San Diego, GmbH)

For more information on our

LymphoTrack Clonality Suite

please contact us at:

[sales@invivoscribe.com](mailto:sales@invivoscribe.com) or

[sales-eu@invivoscribe.com](mailto:sales-eu@invivoscribe.com)

if located outside of North America.

Thank you very much!

Questions?