

# ONCOLOGY: Molecular Testing in NSCLC – Clinical Aspects in Small Specimen Processing

## Pre-Procedural Evaluation

- Choose the best biopsy method to optimize yield (EBUS-TBNA for large mediastinal adenopathy, TTNA for peripheral nodule, etc.)
- Identify reason for biopsy
  - Initial diagnosis
  - Known diagnosis but need additional tissue for molecular testing
- Optimize pre-procedural imaging to maximize procedural yield

## Specimen Collection

- Image guidance to improve sample acquisition
- Utilize ROSE to confirm adequate tissue for testing needs
- Needle gauge (procedure dependent)
- Number of passes
- Operator skill and technique

## Specimen Handling

- Utilizing ROSE to triage specimen
- Collection of specimen within appropriate media (formalin/non-formalin fixatives)
- Perform additional passes for cell block
- Communicate case details with pathology to optimize specimen triage

## Request for molecular studies and follow up on the biopsy results

Initial biopsy reveals adenocarcinoma

PD-L1 immunohistochemistry

Test for actionable mutations (NGS panel testing favored over individual tests)

Initial biopsy reveals adenocarcinoma, but limited tissue remains after diagnostic workup

Communicate presence of limited testing material to the ordering provider and prioritize testing based on discussion

Consider ordering cell-free DNA test (informative, if positive)

Consider repeat biopsy, communicate “molecular priority” protocol for known diagnosis

Patients progressing on initial EGFR TKI

Test for actionable mutations such as T790M, MET amplification, ERBB2/Her-2 amplification

Cell-free DNA test (informative, if positive), otherwise repeat tissue biopsy

Patient progressing after immunotherapy: biopsies remain experimental in this situation.

See online supplement for references and abbreviations: [www.amp.org/PocketGuides](http://www.amp.org/PocketGuides)



# ONCOLOGY: Molecular Testing in NSCLC – Laboratory Aspects in Small Specimen Processing

